

Answer Questionnaire using Case Note Module

1. To enter answers to the questions in the 'Questionnaire', click on the **Open** button in the **Questionnaire** section on the **Case Note** module.

* **Questionnaire**

WG Short Set of Questions on Disability

No Question Answered

Open

2. Select the question and enter the answer. Click on the **Save & Next** button to save the answer and go to the next question.

Answer Questionnaire

WG Short Set of Questions on Disability

Filter

Reset

1. Do you have difficulty seeing, even if wearing glasses?*

[Answer Required]

2. Do you have difficulty hearing, even if wearing hearing aids?*

[Answer Required]

3. Do you have difficulty walking or climbing stairs?*

[Answer Required]

4. Do you have difficulty remembering things?*

[Answer Required]

5. Do you have difficulty (with self-care) such as dressing, grooming, bathing, etc.?*

[Answer Required]

6. Using your usual (customary) language, do you have difficulty understanding or being understood?*

[Answer Required]

1. Do you have difficulty seeing, even if wearing glasses?*

☒ No - no difficulty
 ☐ Yes - some difficulty
 ☐ Yes - a lot of difficulty
 ☐ Cannot do at all

Save

Save & Next

Next

Close Popup

3. After clicking on the **Save** button, a success message will be displayed indicating that the answer has been saved. You may edit the answer or click on the **Delete Answer** button to remove the answer. Click on the **Close Popup** button to go back to the Case Note form.

Answer Questionnaire

WG Short Set of Questions on Disability

Filter

Reset

1. Do you have difficulty seeing, eve...
[Answered]

2. Do you have difficulty hearing, ev...
[Answered]

3. Do you have difficulty walking or c...
[Answered]

4. Do you have difficulty rememberin...
[Answered]

5. Do you have difficulty (with self-c...
[Answered]

6. Using your usual (customary) lang...
[Answered]

Successfully Saved

6. Using your usual (customary) language, do you have difficulty communicating, for example understanding or being understood?*

☒ No - no difficulty
 ☐ Yes - some difficulty
 ☐ Yes - a lot of difficulty
 ☐ Cannot do at all

Prev

Delete


Save

Close Popup

4. After clicking on the **Close Popup** button, you will receive a reminder message to save the Case Note after closing the popup window. Click on the **OK** button to proceed to the Case Note.

secure.therapglobal.net says

Please remember to Save/Update the Case Note form after closing the pop-up.



The saved answers of the questions will appear in the 'Questionnaire' section. If the Questionnaire section is required, then you need to answer the required questions of the Questionnaire before saving/submitting the Case Note.

Questionnaire

WG Short Set of Questions on Disability

Question	Answer
1. Do you have difficulty seeing, even if wearing glasses?	No - no difficulty
2. Do you have difficulty hearing, even if using a hearing aid?	Yes – some difficulty
3. Do you have difficulty walking or climbing steps?	No – no difficulty
4. Do you have difficulty remembering or concentrating?	No – no difficulty
5. Do you have difficulty (with self-care such as) washing all over or dressing?	No – no difficulty
6. Using your usual (customary) language, do you have difficulty communicating, for	No – no difficulty

Open

Note:

Once done answering the questionnaire, don't forget to click on the **Save/Submit** button located at the bottom of the **Case Note** form according to your need.

Cancel

Back

Save

Submit