

Health Tracking Report

Therap's Health Tracking Report module allows users to produce a report for an individual for a chosen time period. The report will contain all the Health related data that has been previously entered for the individual for the selected date range.

Users can generate two types of Health Tracking Reports:

- [Monthly](#)
- [Detailed](#)

Monthly

The Monthly report contains high and low Vital Signs, along with counts of Height/Weight, Blood Glucose, Medications, Seizures, Wounds, Menses, Intake/Elimination and Respiratory Treatments. Each number is then a link to a list of those entries.

1. Click on the **Monthly** link beside Health Tracking Reports under Health Tab.

Health	Data Record Data Due Medications Search Report		SComm ⋮
Agency	<div>Health Tracking</div> <hr/> <div>Appointments New Search Calendar View</div> <hr/> <div>Health Tracking Reports Monthly Detailed</div>		Inbox (4)
Billing			Sent Items
Admin			Compose
			Drafts
			Custom User Group
			Message Audit
			Delete Message

- Users can also click on the **Monthly** link in the Health Tracking Reports dropdown list under Health Tracking, from the **Modules** section of the **Individual Home Page**, by searching an Individual name.

Modules

- ▶ Attendance
- ▶ Case Note
- ▶ General Event Reports (GER)
- ▶ GER Resolution
- ▼ Health Tracking
 - ▶ Appointments
 - ▶ Blood Glucose
 - ▶ Health Care Report
 - ▶ Health Tracking Reports
 - Monthly
 - Detailed
- ▶ Health Tracking Review

2. On the **Health Tracking Monthly Report Page**, Select an **Individual Name** and select **Date range** as required. Then Click on the **Generate** button to see the reports.

Health Tracking Monthly Report

* Individual

Aysha Begum / 987456321



▶ * Begin Date

02/01/2022



▶ * End Date

03/31/2022



Cancel

Generate

- The **Health Tracking Monthly Report page** will appear. Users will be able to see the name of the reports and also the counts of the selected months.
- Users can export the monthly Health Tracking Reports to **PDF** or **excel**.

Health Tracking Monthly Report

Individual: Begum, Aysha
Date: 02/01/2022 - 03/31/2022

	Feb, 22	Mar, 22
Temperature above 101	0	2
Temperature below 97	4	1
Pulse below 60	0	1
Pulse above 100	3	2
Respiration below 12	1	1
Respiration above 20	5	4
Blood pressure high	200/90	150/90
Blood pressure low	110/90	120/80
Weight high	130	130
Weight low	130	130
Bowel Movements	1	2
Laxative	1	1
Suppository	0	1
Enema	0	0
Emesis	1	2
No Intake/Elimination Record	0	0
Menses	0	0
Seizures	5	1
Decubitus	0	0
Other Wound Types	1	0
Number of Medication Started	0	0
Number of Medication Stopped	0	0
Blood Glucose High	6	13
	6	9.5
	6	6

Export to PDF
Export to Excel

[Printables and Exports](#)

Cancel

Back

3. By clicking on the counts, details will be shown in a different window.

Blood pressure high	200/90
Blood pressure low	110/90

Blood pressure low

Individual: Begum, Aysha

Form ID	Blood Pressure (mmHg)	Date	Entered By	Time Zone
HTV-DEMOTPHL-L5X4NJ5XCMUL6	110/90	02/07/2022	Akter, Amina / Academic Coordinator	Asia/Manila

Showing 1 to 1 of 1 entries

4. Users can also see the approved form by clicking on the **Form ID**

Form ID	Blood Pressure (mmHg)	Date	Entered By	Time Zone
HTV-DEMOTPHL-L5X4NJ5XCMUL6	110/90			

Vital Signs Approved

General Information

Individual: Aysha Begum

Program: Training Program

* Reported By: Amina Akter, Academic Coordinator

* Date: 02/07/2022

Notification Level: Low

Vital Signs Information

NOTE: You can view the normal ranges by clicking the "i" icons below

Temperature

Value:

* Time: 07:05 pm

Copy To All

[\[Back To Top\]](#)

Detailed

The detailed Health Tracking report helps users to see any selected form in detail.

1. Click on the **Detailed** link beside Health Tracking Reports under the Health Tab.

Health	Data	Record Data Due Medications Search Report	SComm
Agency	Health Tracking		Inbox (4)
Billing	Appointments	New Search Calendar View	Sent Items
Admin	Health Tracking Reports	Monthly Detailed	Compose
			Drafts
			Custom User Group
			Message Audit
			Delete Message

- Users also can access the **Detailed** link from the 'Modules' section of the **Individual Home Page**, under the **Health Tracking Reports** dropdown list.

Modules

- ▶ Attendance
- ▶ Case Note
- ▶ General Event Reports (GER)
- ▶ GER Resolution
- ▼ Health Tracking
 - ▶ Appointments
 - ▶ Blood Glucose
 - ▶ Health Care Report
 - ▼ Health Tracking Reports
 - Monthly
 - Detailed
- ▶ Health Tracking Review
- ▶ Height/Weight

2. Select the **Individual** name and select the date range as needed. Users can also select the form type(s) from the **Select Form Type(s)** dropdown list. Click on the **Generate Button** to generate the report.

Health Tracking Detailed Report

* Individual

Aysha Begum

Q

→ * Begin Date

02/01/2022

→ * End Date

02/28/2022

→ Select Form Type(s)

3 items selected

Appointments ✓

Blood Glucose ✓

Height/Weight ✓

Immunization

Infection Tracking

Intake/Elimination

Lab Test Result

Medication History

Menses

Respiratory Treatment

Seizures

Skin/Wound

Vital Signs

Cancel

Generate

- The **report** generation date range available for a **Health Care Detailed Report** is 13 months or less.

3. Details of all forms under the health tracking modules, or for the specific form types selected from the dropdown, will be shown on the '**Health Tracking Detailed Report**' page. Users can export the report to a PDF or Excel format through the **Printables and Exports** links.

Health Tracking Detailed Report

Individual: Begum, Aysha
Date: 02/01/2022 - 02/28/2022

Appointments

Form ID	Program	Date	Begin Time	End Time	Appointment With	Specialty	Appointment Status	Address	Phone	Reported By	Time Zone
HTA-DEMOTPHL-L4N4NB7ZJMULU	Training Program (Training Site)	02/19/2022	12:00 pm		Jessia ahmed / Dr.	Dentist	Scheduled	111/J, Road-5, Dhaka-1212, Bangladesh	1212313133	Akter, Amina / Academic Coordinator	Asia/Manila

Blood Glucose

Form ID	Program	Date of Reading	Time	Value	Fasting?	Time Since Last Meal (Hrs)	Insulin Given?	Insulin Amount	Treatment Type	Medication Name	Nurse/Doctor Notified?	Reported By	Time Zone
HTBG-DEMOTPHL-L4Q4NBGWTMULE	Training Program (Training Site)	02/22/2022	03:00 pm	6	No	12 Hrs	Yes	10	Medication		No	Akter, Amina / Academic Coordinator	Asia/Manila
HTBG-DEMOTPHL-L5X4NJ5XAMUL9	Training Program (Training Site)	03/29/2022	06:54 pm	13								Akter, Amina / Academic Coordinator	Asia/Manila

Height/Weight

Form ID	Program	Weight (lb)	Height (Feet)	BMI	Date	Reported By	Time Zone
HTW-DEMOTPHL-L4N4NB8ZX MUL8	Training Program (Training Site)	130	5'	25.39	02/20/2022	Akter, Amina / Academic Coordinator	Asia/Manila

Printables and Exports

4. To open the form, select a row from the list.

Height/Weight

Form ID	Program	Weight (lb)	Height (Feet)	BMI	Date	Reported By	Time Zone
HTW-DEMOTPHL-L4N4NB8ZXMUL8	Training Program (Training Site)	130	5'	25.39	02/20/2022	Akter, Amina / Academic Coordinator	Asia/Manila

Height/Weight Approved

General Information

Individual Aysha Begum

Program Training Program

*** Reported By** Amina Akter, Academic Coordinator

*** Date** 02/20/2022

Notification Level Low

Height/Weight Information

*** Weight (lb)** 130

Height 5' 0"

Feet

Inch

BMI 25.39

[\[Back To Top\]](#)