

## Enter Respiratory Treatment

The Therap Health Tracking **Respiratory Treatment** module allows the users to enter details of Pulse, Respiration Rate, Oxygen Saturation and Lung Sounds before and after a **respiratory treatment**.

Users will need the **HT Submit** role to enter **Respiratory Treatment** information for individuals.

**1. Click on the New link beside the Respiratory Treatment option under the Health tab.**

| To Do          | Health Tracking       |                              |
|----------------|-----------------------|------------------------------|
| Individual     | Appointments          | New   Search   Calendar View |
| Health         | Blood Glucose         | New   Search   Report        |
| Agency         | Height/Weight         | New   Search   Report        |
| Admin          | Respiratory Treatment | New   Search                 |
| Agency Reports | Seizures              | New   Search                 |
|                | Skin/Wound            | New   Search                 |

**2. Select a program from the 'Select Program For Respiratory Treatment' page.**

## Select Program For Respiratory Treatment

All
A
B
C
D
E
F
G
H
I
J
K
L
M
N
O
P
Q
R
S
T
U
V
W
X
Y
Z

Filter

15
Records

| Program Name        | Site Name                                    | Program Type     | Cost Center Number | Program ID |
|---------------------|--|------------------|--------------------|------------|
| Academy             | Therap Global Demonstration Provider - Dhaka | Preschool        |                    | 02         |
| Behavioral Therapy  | DEMO SITE (BD)                               | Behavior Therapy |                    | 07         |
| Bluestar(Day shift) | DEMO SITE (Dhaka)                            | Day Program      |                    | 01         |
| Child Care Center   | Demo-Bashundhara                             | School           |                    | 01         |

3. Select your preferred **Individual** from the Individual List.

## Select Individual For Respiratory Treatment

Program: Behavioral Therapy (DEMO SITE (BD))

All
A
B
C
D
E
F
G
H
I
J
K
L
M
N
O
P
Q
R
S
T
U
V
W
X
Y
Z

Filter

15
Records

| Last Name | First Name | Individual ID | Birth Date |
|-----------|------------|---------------|------------|
| Ahmed     | Tafsir     |               | 02/02/2012 |
| Frances   | Paul       |               | 02/28/2010 |
| Rahman    | Hasan      |               | 04/06/2010 |
| Shonza    | Juliana    |               | 02/02/2010 |
| Uptiti    | Prasant    |               | 03/03/2009 |

**Respiratory Treatment** page will be opened to enter information.


## Respiratory Treatment New ⓘ

### General Information

**Individual** Hasan Rahman

**Program** Behavioral Therapy

**\* Reported By**

**\* Date**  

**Notification Level**

### Respiratory Treatment Information

**\* Event Time**  

#### Pulse

**Before Treatment**

30 - 225

**After Treatment**

30 - 225

#### Respiration Rate

**Before Treatment**

5 - 80

**After Treatment**

5 - 80

#### Oxygen Saturation

**Before Treatment**

95 - 100

**After Treatment**

95 - 100

#### Lung Sounds

**Before Treatment**

**After Treatment**

**Does Order Exist for Chest  
Physiotherapy?** ☐ Yes ☐ No


### Comments

About 3000 characters left

4. On the '**Respiratory Treatment**' page, enter the required information in the 'General Information' section.

## Respiratory Treatment New ⓘ

### General Information

|                    |  |
|--------------------|--|
| Individual         | Hasan Rahman   |
| Program            | Behavioral Therapy   |
| * Reported By      | Imran Rahman, Teacher ▼  |
| * Date             | 05/05/2020  |
| Notification Level | Low ▼  |

5. You can enter information in the Respiratory Treatment Information section. You can change the **Event Time** according to your needs. Information regarding the **Pulse**, **Respiration Rates** and **Oxygen Saturation** can be also entered. The normal ranges for 'Pulse', 'Respiration Rates', 'Oxygen Saturation' are shown under each text field.

### Respiratory Treatment Information

\* Event Time

02:58 pm



#### Pulse

Before Treatment

50

30 - 225

After Treatment

80

30 - 225

#### Respiration Rate

Before Treatment

15

5 - 80

After Treatment

23

5 - 80

#### Oxygen Saturation

Before Treatment

95

95 - 100

After Treatment

102

95 - 100

For the entered values that are out of range, it will show a warning mentioning: *The value entered is not within the normal range for 'Specific' Treatment.* Click on the **Ok** button to continue submitting the form.

The value entered is not within the normal range (95 - 100) for **Oxygen Saturation After Treatment**

Ok

Lung Sounds can be entered if needed.

**Lung Sounds**

|   |                   |
|---|-------------------|
| Before Treatment                          | - Please Select - |
| After Treatment                           |                   |
| Does Order Exist for Chest Physiotherapy? |                   |

**Lung Sounds**

|                            |   |
|----------------------------|---|
| Before Treatment           | Rales   |
| After Treatment            | Clear   |
| Does Order Exist for Chest | <input type="radio"/> Yes <input checked="" type="radio"/> No |

**6.** Click on the **Save** or **Submit** button once you are done entering the required information. If you save the form you will be able to locate it from your Worklist for 30 days.

Cancel

Back

Save

Submit

A confirmation message will show upon successfully submitting the form.

HTR-DEMOTPHL-J773L2BZTMUL6 form has been successfully submitted

Back to [Form](#)