


Search and Update Height and Weight

In order to search and update **Height/Weight** information of an individual, the user will need to be assigned with the **HT View** and **HT Update** role.

1. Click on the **Search** link in the **Height/Weight** row under the **Health** tab.



The screenshot displays the Therap Global interface. On the left, there is a navigation sidebar with a 'Tabs' section containing 'To Do', 'Favorites', 'Individual Home Page', 'Individual', 'Health' (highlighted), 'Admin', and 'Settings'. The main content area is titled 'Health Tracking' and contains a list of health-related categories. The 'Height/Weight' row is highlighted, and the 'Search' link in this row is enclosed in a yellow box. Other categories include Appointments, Blood Glucose, Immunization, Infection Tracking, Intake/Elimination, Lab Test Result, Menses, and Respiratory Treatment, each with 'New' and 'Search' links.

Category	Links
Appointments	New Search Calendar View
Blood Glucose	New Search Report
Height/Weight	New Search Report
Immunization	New Search
Infection Tracking	New Search Report
Intake/Elimination	New Search Report
Lab Test Result	New Search Report
Menses	New Search
Respiratory Treatment	New Search

2. On the **Height/Weight Search** page, enter the necessary parameters to narrow down search results and click on the **Search** button.

Height/Weight Search

Form ID or Individual or Date Range Required

Form ID

Program (Site)

Individual

Entered By

Submit Date From

To

Status

Notification Level

 Clear Selection

Cancel

Search

3. On the Height/Weight Search results page click anywhere on a row of an **Height/Weight** form to open it.

Height/Weight Search

Filter 15 Records

Form ID	Program	Site	Individual	Entered By	Date	Weight (lb)	Height (Feet)	BMI	Status	Notification Level	Time Zone
HTW-DEMOTPHL-HBT4NGYW9MULP	Day Shift	DEMO SITE (BD)	Jones, Steve	Ibrahim, Amit / Supervisor	09/25/2019	130	5' 0"	25.39	Approved	Low	Africa/Kampala

Showing 1 to 1 of 1 entries Previous 1 Next


[New Search](#)
[Excel Exports](#)

[Cancel](#) [Back](#)

4. On the **Height/Weight** page make changes on the form as necessary and click on the **Submit** button to save the changes.

Height/Weight Approved ?

General Information

Individual	Steve Jones
Program	Day Shift
* Reported By	<input type="text" value="Amit Ibrahim, Supervisor"/>
* Date	<input type="text" value="09/25/2019"/> 
Notification Level	<input type="text" value="Low"/>

Height/Weight Information

* Weight (lb)	<input type="text" value="130"/>
Height	<input type="text" value="5"/> <input type="text" value="0"/>
	Feet Inch
BMI	25.39

Comments

Lost five pounds since last measured.

About 2963 characters left

 [Display PDF](#)

- If you would like to send the form via SComm to another user, you may click on the **SComm** button.
- Click the **Delete** button if you need to move the form to a deleted status.