

Search and Update Height and Weight

In order to search and update **Height/Weight** information of an individual, the user will need to be assigned with the **HT View** and **HT Update** role.

1. Click on the **Search** link in the **Height/Weight** row under the **Health** tab.



The screenshot shows the Therap Global interface. On the left, there is a sidebar with a 'Tabs' section containing 'To Do', 'Favorites', 'Individual Home Page', 'Individual', 'Health' (selected), 'Admin', and 'Settings'. The main content area is titled 'Health Tracking' and contains a table of health tracking categories. The 'Height/Weight' row is highlighted, and the 'Search' link in that row is enclosed in a yellow box.

Health Tracking	
Appointments	New Search Calendar View
Blood Glucose	New Search Report
Height/Weight	New Search Report
Immunization	New Search
Infection Tracking	New Search Report
Intake/Elimination	New Search Report
Lab Test Result	New Search Report
Menses	New Search
Respiratory Treatment	New Search

2. On the **Height/Weight Search** page, enter the necessary parameters to narrow down search results and click on the **Search** button.

Height/Weight Search

Form ID or Individual or Date Range Required

Form ID

Program (Site)

Search

Individual

Steve Jones



Entered By

Search

Submit Date From

09/01/2019



To

10/30/2019



Status

- Please Select -



Notification Level

- Please Select -



Clear Selection

Cancel

Search

3. On the Height/Weight Search results page click anywhere on a row of an **Height/Weight** form to open it.


Height/Weight Search

15 Records

Form ID	Program	Site	Individual	Entered By	Date	Weight (lb)	Height (Feet)	BMI	Status	Notification Level	Time Zone
HTW-DEMOTPHL-HBT4NGYW9MULP	Day Shift	DEMO SITE (BD)	Jones, Steve	Ibrahim, Amit / Supervisor	09/25/2019	130	5' 0"	25.39	Approved	Low	Africa/Kampala

Showing 1 to 1 of 1 entries

Previous 1 Next

 New Search
 Excel Exports

4. On the **Height/Weight** page make changes on the form as necessary and click on the **Submit** button to save the changes.


Height/Weight Approved ⓘ

General Information

Individual Steve Jones

Program Day Shift

*** Reported By**

*** Date** 

Notification Level

Height/Weight Information

*** Weight (lb)**

Height
Feet Inch

BMI 25.39

Comments

Lost five pounds since last measured.

About 2963 characters left

 [Display PDF](#)

- If you would like to send the form via SComm to another user, you may click on the **SComm** button.
- Click the **Delete** button if you need to move the form to a deleted status.