

Search MAR Data

Users assigned with the **MAR Data View** role are able to search the MAR data forms.

1. Click on the **Search** link beside the **Data** option under the **Medication Administration Record** section of the **Health** tab.



2. On the **MAR Data Search** page, enter the necessary search parameters and click on the **Search** button.

MAR Data Search

Form ID

→ Individual

Angela Mary



→ Month

5 items selected



February



March



April



May



June



→ Year

2023



2023



Clear Selection

Cancel

Search

3. According to your search criteria, a list of the MAR forms will appear in the search results. To view the MAR form, click on a row.

MAR Data Search

15 Records

Form ID	Individual	Month	Year	Time Zone
MAR-DEMOTPHL-M7Z4NE6Z8833Q	Mary, Angela	June	2023	Asia/Manila
MAR-DEMOTPHL-M7R4NEKYG842G	Mary, Angela	April	2023	Asia/Manila
MAR-DEMOTPHL-M7R4NEKYD842L	Mary, Angela	March	2023	Asia/Manila
MAR-DEMOTPHL-M7K4NEK2284ZD	Mary, Angela	May	2023	Asia/Manila
MAR-DEMOTPHL-M664MNMCP833V	Mary, Angela	February	2023	Asia/Manila

Showing 1 to 5 of 5 entries

< 1 >

New Search

Export to Excel

Cancel

Back

Note: The search results can be exported to an Excel file by clicking on the **Export to Excel** link.

4. Users assigned with the **MAR Data Update** role will be able to update the form by selecting the **Switch to Detail Mode** option on the **Medication Administration Record Data** page. To save it, click on the **Save** button.

Medication Administration Record Data - May , 2023

Individual Name: Angela Mary
Birth Date: 05/01/2010
Created By: Mariam Hossain, Special Educator on 05/17/2023 12:39 PM
Approved By: Mariam Hossain, Special Educator on 05/31/2023 04:46 PM
Form ID: MAR-DEMOTPHL-M7K4NEK2284ZD
Time Zone:Asia/Manila



Legend

■ Administered(New) ■ Administered ■ Missed/ Refused ■ Due
■ LOA(Leave of absence) ■ On hold ■ Deleted ■ User with no Initial

T-Notes

[Add T-Note](#)

[Filter Medications](#) [Record Bulk Data](#) [Open This Form in Configuration Mode](#)

Scheduled Medication(s)

Astemizole - Tablet, Oral (mouth), Scheduled (Medication) [Switch to Detail Mode](#) [Jump to](#)

Strength: 5mg **Prescriber:** Dr. Jhon Paul / Physiotherapist (Demo Hospital Philippines)
Discontinued by: Mariam Hossain on 06/01/2023

Give Amount / Quantity: Once a day Tablet **Frequency:** 1 X DAILY
Begin Date & Time: 11/21/2022 2:00 pm
Schedule Repeat: Every Day 1 time(s) a day **Schedule Time Slot(s):** 2:00 pm
Instruction: Give medicine once a day, 0+1+0

Time	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
	MON	TUE	WED	THU	FRI	SAT	SUN	MON	TUE	WED	THU	FRI	SAT	SUN	MON	TUE	WED	THU	FRI	SAT	SUN	MON	TUE	WED	THU	FRI	SAT	SUN	MON	TUE	WED
2:00 pm	MHN	MHN	MHN	MHN	LOA	LOA	MHN	MHN	MHN	OH	MHN	MHN	MHN	MHN	MHN	OH	MHN	MHN	MHN	MHN	MHN	MHN	MHN	MHN	MHN	MHN	MHN	MHN	MHN	MHN	MHN

Indication/Purpose

To cure, halt, or prevent Allergic rhinitis, unspecified.

Instruction/Comments

Give the medicine once at day.

Fexofenadine - Tablet, Oral (mouth), Scheduled (Medication) [Switch to Detail Mode](#) [Jump to](#)

Strength: 120mg **Prescriber:** Sharma Prendy

Give Amount / Quantity: Once a day Tablet **Frequency:** 1 X DAILY
Begin Date & Time: 04/27/2023
Schedule Repeat: Every Day 1 time(s) a day **Schedule Time Slot(s):** Unspecified
Instruction: Give once daily.

Time	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
	MON	TUE	WED	THU	FRI	SAT	SUN	MON	TUE	WED	THU	FRI	SAT	SUN	MON	TUE	WED	THU	FRI	SAT	SUN	MON	TUE	WED	THU	FRI	SAT	SUN	MON	TUE	WED
	MHN	MHN	MHN	OH	LOA	MHN	MHN	MHN	MHN	MHN	MHN	MHN	MHN	MHN	MHN	LOA	LOA	MHN	MHN	MHN	MHN	MHN	MHN	MHN	MHN	MHN	MHN	MHN	MHN	MHN	MHN

Indication/Purpose

Medication given due to allergy.

Instruction/Comments

Give one tablet everyday.

Other Medication(s)

Neutrogena Sheer Zinc Dry Touch Mineral Face Sunscreen - Creams, Other, Other [Switch to Detail Mode](#) [Jump to](#)

Strength: 50Application **Prescriber:** Dr. Jhon Paul / Physiotherapist (Demo Hospital Philippines)
Give Amount / Quantity: Use whenever going outside or going daylight.Creams **Frequency:** EVERY DAY
Begin Date & Time: 03/08/2023 11:00 am
Schedule Repeat: Every Day 2 time(s) a day **Schedule Time Slot(s):** Unspecified, Unspecified
Instruction: Use whenever going outside or going daylight.

Time	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
	MON	TUE	WED	THU	FRI	SAT	SUN	MON	TUE	WED	THU	FRI	SAT	SUN	MON	TUE	WED	THU	FRI	SAT	SUN	MON	TUE	WED	THU	FRI	SAT	SUN	MON	TUE	WED
	MHN	MHN	MHN	MUR	LOA	MHN	MHN	MHN	MHN	MHN	MHN	OH	MUR	MUR	MHN	MHN	MHN	MHN	OH	MHN	MHN	MHN	LOA	MHN	MHN	MHN	MUR	MUR	MHN	MHN	MHN
	MHN	MHN	MHN	MHN	OH	MHN	MHN	MHN	MHN	MHN	MHN	MHN	MHN	MHN	MHN	MHN	MHN	MHN	LOA	MHN	MHN	MHN	LOA	MHN	MHN	MHN	MUR	MUR	MHN	MHN	MHN

Instruction/Comments

Apply in the outer areas of body whenever going out in daylight.

Allergies

Allergy Status : Known Allergy
Drug Allergy Status : Known Drug Allergy
Allergies : Penicillin

Shared Contact

Name	Organization Name	Specialty	Phone
Dr. Jhon Paul / Physiotherapist	Demo Hospital Philippines	Physical Medicine & Rehabilitation	0112423597

Medical Diagnosis

Diagnosis Coding Type	Diagnosis Code	Description	DSM-5	Billable	Diagnosis Date	Diagnosed By
ICD-10	F84.0 - Autistic disorder	Faces difficulties to communicate and understand what people are saying to her.	Yes	Yes	02/01/2017	Ahmad Imam / Dr. (Demo Hospital)

Dietary Guidelines

Mary can eat only processed foods



[Back](#) [Cancel](#)

[Save](#)

5. Users will receive a success message.

The form MAR-DEMOTPHL-M7Z4NE6Z78334 has been successfully saved