

Hanapin at Mag-Update ng Height at Weight

Upang makapaghanap at mag-update ng impormasyon tungkol sa **Height/Weight**, ang user at dapat nabigyan ng **HT View** at **HT Update** na role.

1. I-click ang link ng **Search** sa hilera ng **Height/Weight** sa **Health** tab.

To Do	Health Tracking	
Individual	Appointments	New Search Calendar View
Health	Blood Glucose	New Search Report
Agency	Height/Weight	New Search Report

2. Sa **Height/Weight** Search page, i-enter ang mga nararapat na datos upang mapaliit ang resulta matapos nito ay i-click ang **Search** button.

Height/Weight Search

Form ID or Individual or Date Range Required

Form ID

Program (Site)

Search

Individual

Steve Jones



Entered By

Search

Submit Date From

09/01/2019



To

10/30/2019



Status

- Please Select -



Notification Level

- Please Select -



Clear Selection

Cancel

Search

3. Sa lumabas na resulta ng **Height/Weight** Search results page, i-click ang anumang nais na hilera upang ito ay buksan.

Height/Weight Search

15 Records

Form ID	Program	Site	Individual	Entered By	Date	Weight (lb)	Height (Feet)	BMI	Status	Notification Level	Time Zone
HTW-DEMOTPHL-HBT4NGYW9MULP	Day Shift	DEMO SITE (BD)	Jones, Steve	Ibrahim, Amit / Supervisor	09/25/2019	130	5' 0"	25.39	Approved	Low	Africa/Kampala

Showing 1 to 1 of 1 entries
Previous 1 Next

[New Search](#)

[Excel Exports](#)

[Cancel](#) [Back](#)

4. Sa page ng **Height/Weight** palitan ang mga nais at i-click ang **Submit** upang ma-save ang mga pinalitan.

Height/Weight Approved ⓘ

General Information

Individual	Steve Jones
Program	Day Shift
* Reported By	Amit Ibrahim, Supervisor
* Date	09/25/2019
Notification Level	Low

Height/Weight Information

* Weight (lb)	130	
Height	5	0
	Feet	Inch
BMI	25.39	

Comments

Lost five pounds since last measured.

About 2963 characters left

 Display PDF

Cancel Back

SComm Delete Submit

- Kung nais na ipadala ang form na ito gamit ang Scomm sa ibang user, i-click lamang ang **Scomm** button.
- I-Click ang **Delete** button kung nais i-delete ang form na ito.