

Mag-Generate ng ISP Report - Programmatic Report

Mag-generate ng ISP Reports sa mga serbisyong ibinigay sa mga indibidwal upang masundan ang pagbigay ng serbisyo at ang progress ng kanilang goals.

Ang mga users na may role na **ISP Report** ay makakagawa ng report ng ISP Program sa mga indibidwal na nasa kanilang caseload. Makakagawa ng reports sa mga approved at discontinued ISP Programs.

Program: Profile: Module:	No Program Selected Initial Search	Choose Program
To Do	Care	Issue Tracking Wy Issues
Individual	T-Log New Search Archive	SComm "
Health	Case Note New Search Archive Bulk PDF	Inbox Sent Items
Agency	ISP Data New Search Report Search Report Archive	Compose Drafts Custom User Group
Agency	ISP New Search Review Acknowledge Archive	Classes
Reports	ISP Program New Search Archive Acknowledgement Report	Overdue Due
Home Page	ISP Program Template Library New Draft Published Search	View Sign ups View Results/Notes

1. I-click ang **Report** link na makikita sa hilera ng **ISP Data** sa Individual Tab.

2. Piliin ang program mula sa Select Program sa page ng ISP Program.



Select Program For ISP Program

AII A B C D E F G	H I J K L M N O P Q	R S T U V	W X Y Z	
Filter			15	~ Records
► Program Name	\$	Program Type	Cost Center \$\\$ Number	Program
Child Care Center	Care Home	School		01
Day Shift	Special School	School		01
Demo Program	Demo Site	School		555
Chausing 1 to 2 of 2 antrias				

3. Piliin ang indibidwal mula sa Individual List sa page ng ISP Program.

Individual List for ISP Program	
---------------------------------	--

AII A B C D E F G	H I	J K L	M N	0 P	Q R	S	T U	۷V	v x	Y	Ζ	
Filter											15	 Records
Last Name	•	First Name			ŧ	•	Individu	al ID				¢
Tareq		Taufiq										

4. Piliin ang ISP Program name at i-click ang **Next** button upang makagenerate ng report.

Person-Centered. Data-Driven. **ISP Program List** Filter Records You have selected 1 items. 15 Program Site **Program Creation** Approve Status Individual ISP Program Form ID Name Name Entered By Date Date Approved **•** Time Zone ISP-DEMOTPHL-Demo Sonia Ferdaus, Therap 05/01/2019 07/22/2019 Asia/Dhaka 1 Demo Tareq, Speech Therapy Approved H9Q4NFNYTMULU Program Site Taufiq Plan Admin Showing 1 to 1 of 1 entries 1 Next Previous Cancel

nerap Global

Mayroong apat na uri : **Programmatic Report, Clinician Report, Data Collection Monthly Report, at Hab Documentation Record.**

5. Naka display ang page ng Report Criteria. Piliin ang nais sa field at i-click ang Generate button. Upang maka-generate ng Programmatic report, piliin ang Programmatic mula sa Report Type field.

Report Criteria	
Report Type	Programmatic ~
Format	Table ~
Time Zone	Asia/Dhaka
Start Date	06/01/2019
End Date	06/30/2019
Review Period	Month
Include Linked ISP Program(s)	⊖ Yes ● No
Show Deviation from Baseline Score	○ Yes
Show Deviation from Previous Score	● Yes ◯ No
Score Showing Option	Both
Add progress and comment per ISP Program	○ Yes ● No

Terap Global

Centered. Data-Driven.

Selected ISP Program			
Form ID	ISP Program	Approve Date	
ISP-DEMOTPHL-H9Q4NFNYTMULU	Speech Therapy Plan	07/22/2019	Select Task/Score
Back			Generate

Naglalaman ang page ng **Report Criteria** ng mga sumusunod na fields:

Report Type: Programmatic, Clinician, Data Collection Monthly, o Hab Documentation Record. Upang mag-generate ng Programmatic Report piliin ang Programmatic Report mula sa drop down menu.

Format: Maaaring mamili mula sa Table, Bar Graph, o Line Graph para sa report format. (Ito ay available lamang kapag ang napili ay Programmatic report type.)



Start Date & **End Date**: Pumili ng date range ng report. Sa Programmatic Reports, ang date range ay may maximum na 2 taon.

Review period: Mamili kung paano gusto ipakita ang mga datos pwedeng ayon sa Day, Week, Month, Quarter, Year, o pwede rin ang Custom option.

Score Showing Option: Maaaring mamili ng Percentage only, Trial Count only, o Both. (Ito ay available lamang kapag ang napili ay Programmatic report type.)

Programmatic Report

Ito ang mga nilalaman ng **Programmatic Report**.

- Sa Report Description field, maaaring maglagay ng short name para sa report, ito ay kailangan kung nais mong i-save ang report.
- Ang individual name, date range, review period at time zone.
- Ang Generated Report section ay magpapakita ng mga datos para sa ISP na nakolekta mula sa napiling date range para sa ISP Program.

Ang mga sumusunod ay ginamit para sa Programmatic Reports:

- **S**: Score ng bawat task.
- $\circ~{\bf C}:$ Trial Count, o kung ilang beses binilang ang task
- **DB**: Deviation sa pagitan ng task score at ng Baseline Score. (Hindi ito kailangan gamitin kung hindi naglagay ng Baseline score)

Person-Centered, Data-Driven.

				Programmatic Report
				Status: New
Report Type : Prog	rammatic			
Report Description *	Taufiq June Repo	ort		
Individual Name	Taufiq Tareq,		Review Period	Month
Date Range	06/01/2019 - 06	/20/2019	Time Zone	Asia/Dhaka
- Generated Report				
			S = Score	e, C = Trial Count, DB = Deviation from Baseline Score, DP = Deviation from Previous Score
Speech Therapy Pla	an			
Form ID		ISP-DEMOTPHL-H9Q4NFNYTMULU		
Status		Approved (Dynamic)		
ISP Program		Speech Therapy Plan		
Goal/Service		Ability to comply with verbal direction	and develop independence	
Criteria for Completion				
Scoring Method		Level of Independence		
Score(s)				
DECL		Declined		
FULL		Full Physical Prompt		
PART		Partial Physical Prompt		
DEMO		Demonstration		
VRBL		Verbal Cue		
INDP		Independence		
Task(s)				
1		Give me Ball		
2		Clap Hands		
3		Say Hello (maintains eye contact for 5	5 sec)	

Mayroon tatlong uri ng format ang **Programmatic Reports** na makikita.

Ito ay ang mga sumusunod : **Table, Bar Graph** at **Line Graph**.

Table

Teek	Desilies	06/01/2019 - 06/30/2019			07/01/2019 - 07/20/2019				Total	
Task	Daseine	S	С	DB	s	с	DB	DP	s	с
1	DECL 0 %	12.5 %	1	12.5 %	0 %	0	0 %	-12.5 %	12.5 %	1
	FULL 0 %	0 %	0	0 %	0 %	0	0 %	0 %	0 %	0
	PART 0 %	0 %	0	0 %	0 %	0	0 %	0 %	0 %	0
	DEMO 0 %	12.5 %	1	12.5 %	0 %	0	0 %	-12.5 %	12.5 %	1
	VRBL 0 %	75 %	6	75 %	0 %	0	0 %	-75 %	75 %	6
	INDP 0 %	0 %	0	0 %	0 %	0	0 %	0 %	0 %	0
	DECL 0 %	12.5 %	1	12.5 %	0 %	0	0 %	-12.5 %	12.5 %	1
2	FULL 0 %	0 %	0	0 %	0 %	0	0 %	0 %	0 %	0
	PART 0 %	0 %	0	0 %	0 %	0	0 %	0 %	0 %	0
	DEMO 0 %	37.5 %	3	37.5 %	0 %	0	0 %	-37.5 %	37.5 %	3
	VRBL 0 %	37.5 %	3	37.5 %	0 %	0	0 %	-37.5 %	37.5 %	3
	INDP 0 %	12.5 %	1	12.5 %	0 %	0	0 %	-12.5 %	12.5 %	1
	DECL 0 %	50 %	4	50 %	0 %	0	0 %	-50 %	50 %	4
	FULL 0 %	25 %	2	25 %	0 %	0	0 %	-25 %	25 %	2
3	PART 0 %	0 %	0	0 %	0 %	0	0 %	0 %	0 %	0
·	DEMO 0 %	12.5 %	1	12.5 %	0 %	0	0 %	-12.5 %	12.5 %	1
	VRBL 0 %	12.5 %	1	12.5 %	0 %	0	0 %	-12.5 %	12.5 %	1
	INDP 0 %	0 %	0	0 %	0 %	0	0 %	0 %	0 %	0
	DECL 0 %	25 %	6	25 %	0 %	0	0 %	-25 %	25 %	6
	FULL 0 %	8.33 %	2	8.33 %	0 %	0	0 %	-8.33 %	8.33 %	2
Tetal	PART 0 %	0 %	0	0 %	0 %	0	0 %	0 %	0 %	0
10101	DEMO 0 %	20.83 %	5	20.83 %	0 %	0	0 %	-20.83 %	20.83 %	5
	VRBL 0 %	41.67 %	10	41.67 %	0 %	0	0 %	-41.67 %	41.67 %	10
	INDP 0 %	4.17 %	1	4.17 %	0 %	0	0 %	-4.17 %	4.17 %	1



Bar Graph





Line Graph





Kung nais ng karagdagang kaalaman tungkol sa Clinician Report i-click dito.
Kung nais ng karagdagang kaalaman tungkol sa Hab Documentation Record i-click dito.
Kung nais ng karagdagang kaalaman tungkol sa Data Collection Monthly i-click dito.